



**CUPID**  
LIMITED

# CupiKIT

**DENGUE IgG /IgM Test Kit**  
For In Vitro Diagnostic Use Only

## PRODUCT SPECIFICATION

- Chromatographic immunoassay
- Detect Dengue Antibodies
- Result interpretation time: 5 minute
- Specimen: Serum/Plasma
- Sensitivity: 100%
- Specificity: 100%
- Storage Temperature : 2°C -40°C

## MATERIALS PROVIDED


### PROVIDED

- Test Device
- Dropper
- Desiccant
- Instructions for Use

### REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

**PACK SIZE: 50 Test Kit**



**Manufactured by: Cupid Limited, A-68, M.I.D.C. (Malegaon),  
Sinnar, Nashik - 422 113, Maharashtra, India,  
Mfg. Lic. No.: MFG/IVD/2021/000122**