



CUPID
LIMITED

CupiKIT

LH OVULATION Test Kit

For In Vitro Diagnostic Use Only

PRODUCT SPECIFICATION

- Chromatographic immunoassay
- Detection of Luteinizing hormone
- Result interpretation time: 5-10 minutes
- Specimen: Urine
- Sensitivity: 99.8 %
- Specificity: 100%
- Storage Temperature : 2°C -40°C

MATERIALS PROVIDED

PROVIDED

- Test Device
- Dropper
- Desiccant
- Instructions for Use

REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer



Manufactured by: Cupid Limited, A-68, M.I.D.C. (Malegaon),
Sinnar, Nashik - 422 113, Maharashtra, India,
Mfg. Lic. No.: MFG/IVD/2021/000122