



CUPID
LIMITED



ACCURACY



EASY



**GET RESULT
IN 5 MINUTES**



PREGNANCY HCG TEST KIT
For In Vitro Diagnostic Use Only

PRODUCT SPECIFICATION

- HCG Test Device for detection of hCG in Urine.
- Easy to use
- High Specificity
- Specimen: Urine
- Sensitivity: 25 mIU/ml
- Storage temperature : 2°C -40°C

MATERIALS PROVIDED

PROVIDED

- Test Device
- Dropper
- Desiccant
- Instructions for Use

REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

PACK SIZE: 100 Test Kit



**Manufactured by: Cupid Limited, A-68, M.I.D.C. (Malegaon),
Sinnar, Nashik - 422 113, Maharashtra, India,
Mfg. Lic. No.: MFG/IVD/2021/000122**