



CUPID
LIMITED

CupiKIT

Typhoid IgG/IgM Test Kit

For In Vitro Diagnostic Use Only

PRODUCT SPECIFICATION

- Chromatographic immunoassay
- Detection of *S. typhi* specific antibodies
- Result interpretation time: 10-15 minutes
- Specimen: Serum / Plasma
- Sensitivity: 100 %
- Specificity: 100%
- Storage Temperature : 2°C -40°C

MATERIALS PROVIDED

PROVIDED

- Test Device
- Dropper
- Desiccant
- Assay Buffer
- Instructions for Use

PACK SIZE: 50 Test Kit



**Manufactured by: Cupid Limited, A-68, M.I.D.C. (Malegaon),
Sinnar, Nashik - 422 113, Maharashtra, India,
Mfg. Lic. No.: MFG/IVD/2021/000122**