

Viral Transport Medium

PRODUCT SPECIFICATION

Intended use: for the collection and transport of clinical specimens containing viruses

- Easy to handle & use
- No other accessories required
- Storage Temperature : 2°C -40°C

KIT CONTENT

DESCRIPTION	QUANTITY	STORAGE
CUPIKIT_ Viral Transport Medium	3ML	Room Temperature
Sterile nylon Swab with breakpoint (Thick for Throat and Thin for Nose sampling)	2	Room Temperature
Instructions for Use	1	—



PACK SIZE 100 Test Kit

