



CupiKIT

DENGUE Ns1 Antigen Test Kit For In Vitro Diagnostic Use Only

PRODUCT SPECIFICATION

- Chromatographic immunoassay
- Detection of Dengue NS1 specific antigens
- Specimen: Serum/Plasma
- Result interpretation time: 15 – 20 mins
- Sensitivity: 100%
- Specificity: 100%
- Storage Temperature : 2°C -40°C

MATERIALS PROVIDED

PROVIDED

- Test Device
- Dropper
- Desiccant
- Instructions for Use

PACK SIZE: 50 Test Kit



Manufactured by: **Cupid Limited**, A-68, M.I.D.C. (Malegaon),
Sinnar, Nashik - 422 113, Maharashtra, India,