



CUPID
LIMITED

CupiKIT

Dengue Combo Rapid Test Kit (Dengue Ns1 Ag & IgG/IgM)

For In Vitro Diagnostic Use Only

PRODUCT SPECIFICATION



- Chromatographic immunoassay
- Detection of Dengue Ns1 Ag & Dengue IgG/ IgM Test Kit.
- Specimen: Serum/Plasma
- Result interpretation time: 15 – 20 mins
- Sensitivity: 100%
- Specificity: 100%
- Storage Temperature : 2°C -40°C

MATERIALS PROVIDED

PROVIDED

- Test Device
- Dropper
- Desiccant
- Instructions for Use

PACK SIZE: 10 Test Kit



**Manufactured by: Cupid Limited, A-68, M.I.D.C. (Malegaon),
Sinnar, Nashik - 422 113, Maharashtra, India,**