



**CUPID**  
LIMITED

# CupiKIT

## HIV / Syphilis Combo Test Kit

For In Vitro Diagnostic Use Only

### PRODUCT SPECIFICATION

- Chromatographic immunoassay
- Detection of HIV and Syphilis antibodies
- Result interpretation time: 20 minutes
- Specimen: Whole blood/ Serum/ Plasma


For HIV:	For Syphilis:
■ Sensitivity - 100 %	■ Sensitivity - 95 %
■ Specificity - 100 %	■ Specificity - 100 %

- Storage Temperature : 2°C -40°C

### MATERIALS PROVIDED

#### PROVIDED

- Test Device
- Dropper
- Desiccant
- Instructions for Use



PACK SIZE: 50 Test Kit



**Manufactured by: Cupid Limited, A-68, M.I.D.C. (Malegaon),  
Sinnar, Nashik - 422 113, Maharashtra, India,**